


LYMPHOEDEMA CATALOGUE 2015



CONTENTS



Pumps & Equipment	2
Pump Garments	5
Compression Garments	11



PUMPS & EQUIPMENT

LYMPHA PRESS® MINI

Treating Lymphoedema correctly means a lifelong commitment to a regimen of care. Home therapy with the Lympha Press® Mini is easy and effective, saves clinic visits, and increases patient compliance for better treatment results.

INDICATIONS

- Lymphoedema of legs and arms.
- Post-mastectomy arm Lymphoedema.
- Healing of venous stasis ulcers.
- Promotion of wound recovery.
- Relief from muscular and nervous pain.
- Relief of discomfort from leg fatigue.

Technical Specifications	
Model	201
Cycle	Sequential
Cycle Time	30 sec
Pressure Range	20-80 mmHg
Size	14.5/38.5/30.5cm
Weight	5.1kg
Electrical	230V / 115V 50/60 Hz

Garment choice and sizes - see separate garment catalogue.

FEATURES FOR LYMPHA PRESS® MINI

- Up to 24 overlapping pressure cells.
- Fully variable pressure range 20 to 80 mmHg for safe use at home.
- User friendly - easy to use.
- Quiet.
- Compact, easy to store.
- Sturdy and reliable.
- Portable.
- Can treat up to two limbs simultaneously.
- Compatible with comprehensive range of quality garments.
- Easy reading controls (Pressure and Treatment duration).

Code

L1000024A



LYMPHA PRESS® OPTIMAL

Based on the reliable and durable Lympa Press® MINI line, the Lympa Press® Optimal is a pressure controlled device and features Pretherapy treatment also for home care patients.

The Pretherapy Treatment is a gentle massaging of the proximal lymphatics, with the purpose of easing the lymph liquids transit, based on MLD (manual lymph drainage) procedures.

The Pretherapy starts the treatment at the proximal end of the treated limb.

Lympa Press® Optimal then moves on to the main treatment, either sequential buildup pressure or Wave Peristaltic massage.

Technical Specifications

Model	1201
Cycle	Sequential Multi Cycles
Pressure Range	20-90mmHg
Size	14.5x38.5x30.5cm
Weight	6.1kg
Electrical	230V / 115V 50/60 Hz

Garment choice and sizes - see separate garment catalogue.

NEW FEATURES

- Pressure controlled device.
- Pretherapy cycle preceding all treatments.
- Pretherapy cycle override option.
- Remote control pause switch.
- Informative LCD screen.
- Quick-connect connector for hose bundles.

OTHER FEATURES

- Up to 24 overlapping air cells.
- Fully variable pressure range 20 to 90 mmHg.
- User friendly - easy to use.
- Quiet.
- Compact, easy to store.
- Sturdy and reliable.
- Portable.
- Can treat up to two limbs simultaneously.
- Compatible with most Lympa Press garments.
- Easy reading controls.

Code

L10002C7



PUMP GARMENTS



COMFY® LEG SLEEVE

Comfy Leg Sleeve is a multiple air-cell compression garment for the treatment of full legs - foot, calf & thigh. Use comfy sleeve with Lympha Press Optimal and Lympha Press Mini.

Code	Size	A	B	C	D	E	F	L
L41L2285	2-850mm	700mm	600mm	470mm	420mm	480mm	260mm	840mm
L41L2295	2-950mm	780mm	680mm	510mm	430mm	480mm	260mm	950mm
L41L2365	3-650mm	810mm	730mm	600mm	560mm	540mm	370mm	650mm
L41L2375	3-750mm	850mm	760mm	610mm	560mm	540mm	370mm	750mm
L41L2385	3-850mm	900mm	800mm	630mm	560mm	540mm	370mm	850mm


1. General tolerance: ± 2 millimeters. 2. A, B, C, D & E dimensions are circumferences.



PUMP GARMENTS


COMFY® LEG SLEEVE PARALLEL EXPANDER

Comfy Leg Sleeve Parallel Expander allows an increase in the circumference of the Comfy Leg Sleeve. Each expander adds 100mm from the top to the bottom of the sleeve.



COMFY® LEG SLEEVE TAPERED EXPANDER

Comfy Leg Sleeve Tapered Expander allows an increase in the circumference of the Comfy Leg Sleeve. Each expander adds 100mm above the knee tapered to 40mm at the ankle.



Code	Use with Sleeve
L415A097	L41L2385
L415A087	L41L2285, L412375
L415A101	L41L2295
L415A078	L41L2265

Code	Size
L415T087	L41L2285, L41L2375
L415T094	L41L2295
L415T078	L41L2365
L415T097	L41L385

LYMPHA® PANTS

Lympha Pants® is designed to give best lymph drainage for both legs and lower abdomen.

Lympha Pants® is also very efficient in cases where only one leg suffers from Lymphoedema.

Use Lympha Pants® with Mego Afek Devices: Lympha Press® Optimal, Lympha Press Mini®.

INSTRUCTIONS FOR WEARING THE LYMPHA PANTS®

1. Arrange the Lympha Pants® on the surface where treatment will take place.
2. Open the main zipper 2/3 & lie down inside the Lympha Pants®.
3. There are three zippers for three different pant sizes. Close the zipper which adapts the size to meet your needs, closing the pants.
4. If the pants are still too small, add Expanders® to adapt the pants size to your size.
5. On each leg, near the ankle, there are three (3) small zippers. Open the one that is closest to the patient's heel. It will allow the foot to remain as close as possible to its normal position (90° to leg) during treatment.

RECOMMENDATIONS

- Use Lympha Pants® while lying down
- Wear cotton training pants under the Lympha Pants

CONNECTING THE HOSE BUNDLE

Insert the hose bundle connectors (in direction pointed by red arrows) into the corresponding air socket receptacles of the device.


DISCONNECTING THE HOSE BUNDLE

Press down the snap tooth to release the hose bundle and gently pull out the bundle from the device.

Code

L4MP2126

LYMPHA PANTS SIZING



PUMP GARMENTS

LYMPHA® PANTS EXPANDER

Lympha Pants Expanders® allow extending the circumference of Lympha Pants®. Use Lympha Pants Expanders® with Mego Afek: Lympha Pants®.

You can use up to three pants expanders on each pant leg. Each expander (right & left) adds about 10cm to the circumference of the pant (leg & abdomen).

Pants Expanders come in sets of two. One for each side of the pants.

INSTRUCTIONS FOR USE

1. Open the zipper on the pant leg that you want to expand.
2. Lay out the Pant Expander® next to one side of the open zipper.
3. Close one side of the Pant Expander® zipper with one side of the pants zipper.
4. Close 1/3 of the other side of the Pant Expander® zipper with the other side of the pants zipper.
5. Use the expanded Lympha Pants® in the standard manner.



Code


L4Z5LP00

Unit

Pair

LYMPHA® JACKET

Designed to give best lymph drainage for both arms and upper abdomen. Very useful and efficient where both or only one arm suffers from Lymphoedema.



Code	Style	Arm 1	Arm 2	Chest Circ.
L4MJ2125	Short	800mm	570mm	1050-1230mm
L4MJ2135	Long	960mm	720mm	1050-1230mm

COMFY® ARM SLEEVE

The Comfy Arm Sleeve 1 - 75 offers a new and greatly improved way to treat Lymphoedema in the upper limbs. The treated area does not stop at the upper arm but now includes the armpit as well as the upper shoulder, one side of the back, chest and abdomen.



Code	Size	Orientation	Torso Circ.
L41LS175L	Small	Universal	1050-1150mm
L41LM175L	Medium	Universal	1110-1250mm
L41LL175L	Large	Universal	1220-1400mm

COMPRESSION GARMENTS



OAPL STOCKING MEASURING GUIDE



Knee Length



Thigh Length

SIZING THE STOCKINGS

Care should be taken in determining the correct size of the stocking.

MEASURING

Knee Length Stockings

- Start by measuring the circumference of the leg at the mid calf 'A'
- Next measure the length of the leg 'H' to 'K' measuring from the heel to the back of the knee. (popliteal fossa)
- Note: Knee Length stockings should finish 3cm below the popliteal fossa

THIGH LENGTH & THIGH WITH WAIST BELT STOCKINGS

- Measure circumference of the leg at the mid calf 'A'
- Next measure the length of the leg from the heel, 'H' to the gluteal fold 'G'
- Note: Thigh Length stockings should end at Gluteal fold.

IMPORTANT

For best results, instructions for applying the stockings should be carefully followed. Failure to do so could result in restricted circulation or reduced effect.

OAPL stock the full range of Jobst Stockings. For a comprehensive Jobst guide please contact customer service or your local representative.

OAPL ANTI-EMBOLISM STOCKINGS

OAPL Anti Embolism stockings are designed to support and aid in venous return, recommended for use by every recumbent patient to prevent DVT Prophylaxis. OAPL Anti Embolism stockings are latex free graduated compression, 18mmHg at the ankle, comprising of 76% Nylon, 24% Lycra Spandex. OAPL Anti Embolism stockings should be used with patients with a history of phlebitis or thrombosis or immobile, aged patients suffering venous insufficiency. Also used post-operatively for oedema management. Stockings can be washed or sterilised (Autoclave at 121 degrees) or hand wash in warm water with mild soap and drip dried.

Code	Size	Variation	Height	Calf
41031J	XS	Regular	Up to 400mm	Under 250mm
41031A	S	Regular	Up to 400mm	250-300mm
41031C	M	Regular	Up to 420mm	300-370mm
41031E	L	Regular	Up to 450mm	370-440mm
41031G	X	Regular	Up to 480mm	440-500mm
41031B	S	Long	Over 400mm	250-300mm
41031D	M	Long	Over 420mm	300-370mm
41031F	L	Long	Over 450mm	370-440mm
41031H	XL	Long	Over 480mm	440-500mm



OAPL ANTI-EMBOLISM STOCKINGS THIGH LENGTH

Code	Size	Variation	Height Gluteal Fold	Calf
41028A	S	Regular	Up to 820mm	250-300mm
41028C	M	Regular	Up to 820mm	300-370mm
41028E	L	Regular	Up to 820mm	370-440mm
41028G	X	Regular	Up to 820mm	440-500mm
41028B	S	Long	Over 820mm	250-300mm
41028D	M	Long	Over 820mm	300-370mm
41028F	L	Long	Over 820mm	370-440mm
41028H	X	Long	Over 820mm	440-500mm



OAPL ANTI-EMBOLISM STOCKING THIGH LENGTH

Code	Size	Variation	Height Gluteal Fold	Calf
41027A	S	Regular	Up to 820mm	250-300mm
41027C	M	Regular	Up to 820mm	300-370mm
41027E	L	Regular	Up to 820mm	370-440mm
41027G	X	Regular	Up to 820mm	440-500mm
41027B	S	Long	Over 820mm	250-300mm
41027D	M	Long	Over 820mm	300-370mm
41027F	L	Long	Over 820mm	370-440mm
41027H	XL	Long	Over 820mm	440-500mm



COMPRESSION GARMENTS

OAPL TRAVEL SOCKS

OAPL Travel Socks minimise the risk of DVT forming in the legs during travel by increasing circulation. The latex free, hospital grade stockings feature graduated compression, 14-18mmHg. Comprising of 76% Nylon, 24% Spandex they are available in Black or Beige, knee length. They are suitable for tired, achy legs (heaviness), ankle swelling and moderate to pronounced spider and varicose veins. Hand wash in warm water with mild soap.

Code	Size	Colour	Height	Calf
41032	Small	Black	Up to 400mm	250-300mm
41033	Medium	Black	Up to 420mm	300-370mm
41034	Large	Black	Up to 450mm	370-440mm
41035	XLarge	Black	Up to 480mm	440-500mm
41036	Small	Beige	Up to 400mm	250-300mm



THERAFIRM MEN'S SOCKS 15-20mmHg

Designed to help relieve minor ankle or leg swelling, minor varicosities and leg fatigue. Soft comfort band is non-restrictive while staying in place all day. Reinforced toe and heel for durability. Sheer, two-way stretch nylon and spandex. Available in closed toe styles. Colors available are sand and black.

Code	Size	Variation	Ankle	Calf	Length
68310	S	Black	180-230mm	270-380mm	280-360mm
68311	M	Black	230-280mm	280-420mm	310-380mm
68312	L	Black	250-360mm	310-430mm	330-410mm
68340	XL	Black	280-380mm	330-480mm	360-460mm
68353	S	Brown	180-230mm	270-380mm	280-360mm
68313	M	Brown	230-280mm	280-420mm	310-380mm
68314	L	Brown	250-360mm	310-430mm	330-410mm
68342	XL	Brown	280-380mm	330-480mm	360-460mm



THERAFIRM WOMEN'S KNEE HIGH 15-20mmHg

Designed to help relieve minor ankle or leg swelling, minor varicosities and leg fatigue. Soft comfort band is non-restrictive while staying in place all day. Reinforced toe and heel for durability. Sheer, two-way stretch nylon and spandex. Available in closed toe styles. Colors available are sand and black.

Code	Size	Variation	Ankle	Calf	Length
68128	S	Black	160-180mm	280-340mm	310-460mm
68124	M	Black	190-220mm	330-390mm	330-480mm
68125	L	Black	240-270mm	380-440mm	360-510mm
68132	XL	Black	280-290mm	450-470mm	380-560mm
68200	S	Sand	160-180mm	280-340mm	310-460mm
68201	M	Sand	190-220mm	330-390mm	330-480mm
68202	L	Sand	240-270mm	380-440mm	360-510mm
68203	XL	Sand	280-290mm	450-470mm	380-560mm



THERAFIRM THIGH HIGH WITH LACE 15-20mmHg

Designed to help relieve minor ankle or leg swelling, minor varicosities and leg fatigue. Over-the-knee style stockings feature non-restrictive lace band with silicone border to keep the stockings in place. Reinforced toe and heel for durability. Made from sheer, two way stretch nylon and spandex. Available in closed toe style.

Code	Size	Variation	Ankle	Calf	Length
68450	S	Sand	160-180mm	280-340mm	580-640mm
68451	M	Sand	190-230mm	410-440mm	640-710mm
68452	L	Sand	230-270mm	440-480mm	710-790mm
68453	XL	Sand	280-290mm	480-530mm	790-860mm
68454	XXL	Sand	300-320mm	530-580mm	860-910mm



THERAFIRM KNEE HIGH 15-20mmHg - CLASS 1

Available in Closed & Open toe. Class 1 stockings are designed to offer moderate to firm support. These stockings may be prescribed for relief from chronic ankle swelling, prophylaxis or clotting after uncomplicated vein stripping, mild Lymphoedema and more.

Code	Size	Variation	Ankle	Calf	Length
79320	S	Open-Beige	160-180mm	270-330mm	280-380mm
79321	M	Open-Beige	190-220mm	290-360mm	310-430mm
79322	L	Open-Beige	220-250mm	320-380mm	330-460mm
79323	XL	Open-Beige	250-290mm	340-410mm	360-480mm
79324	XXL	Open-Beige	230-320mm	370-430mm	410-530mm
79325	XXXL	Open-Beige	320-350mm	420-520mm	410-530mm
79326	XXXXL	Open-Beige	350-380mm	510-610mm	410-530mm
75112	S	Closed-Beige	160-180mm	270-330mm	280-380mm
75113	M	Closed-Beige	190-220mm	290-360mm	310-430mm
75114	L	Closed-Beige	220-250mm	320-380mm	330-460mm
75115	XL	Closed-Beige	250-290mm	340-410mm	360-480mm
75116	XXL	Closed-Beige	230-320mm	370-430mm	410-530mm



THERAFIRM KNEE HIGH 15-20mmHg - CLASS 1

Available in Closed & Open toe. Class 2 stockings are designed to offer extra firm support. These stockings may be prescribed in cases of severe or chronic venous insufficiency related to acute leg swelling, Lymphoedema etc.

Code	Size	Variation	Ankle	Calf	Length
76505	S	Open-Beige	160-180mm	270-330mm	280-380mm
76506	M	Open-Beige	190-220mm	290-360mm	310-430mm
76507	L	Open-Beige	220-250mm	320-380mm	330-460mm
76508	XL	Open-Beige	250-290mm	340-410mm	360-480mm
76509	XXL	Open-Beige	230-320mm	370-430mm	410-530mm
76600K	S	Closed-Beige	160-180mm	270-330mm	280-380mm
76601K	M	Closed-Beige	190-220mm	290-360mm	310-430mm
76602K	L	Closed-Beige	220-250mm	320-380mm	330-460mm



COMPRESSION GARMENTS

OAPL with its long heritage in both circular knit and custom compression garments is proud to act as the exclusive distributor of the Ofa Bamberg range. Ofa Bamberg stands for the most modern technology in combination with more than 80 years of experience in medical compression. They are also proud of the fact that today their product is still exclusively designed and manufactured in Germany.

The quality of our compression products is the first priority. We see a high level of quality assurance as our foremost duty and guarantee to provide, in accordance with the listed certifications, safe and medically correct care.

ofa bamberg

we care for you



AVAILABLE
IN OFF
THE SHELF
& MADE TO
MEASURE

CIRCULAR KNIT

MEMORY

Transparent full-fashioned look. Delicate sheen, velvety soft texture, smooth knitted fabric. Easier to don due to specially developed compression threads. The most modern functional yarns providing for improved comfort and providing optimal fit especially in areas such as the hollow of the knee and at the ankle. Soft, elastic toe and heel. Standard and made-to-measure options.



MEMORY WITH ALOE VERA

Micro encapsulated combination of ingredients – integrated into the yarn. Moisturises through Aloe Vera extracts and Vitamin E. Protects, soothes and cares for the skin. Transparent full-fashioned look. Delicate sheen, velvety soft texture, smooth knitted fabric. Easier to don due to specially developed compression threads. The most modern functional yarns for improved comfort and optimal fit especially in areas such as the hollow of the knee and at the ankle. Soft, elastic toe and heel. Standard and made-to-measure options.



MAX MEDICAL STOCKINGS

Sock style stocking for men. Skin friendly yarn with 50% cotton content. Thermo and moisture regulating. Practical 3-size system. Standard product only.



LASTOFA

Soft, opaque knitted fabric. Fine and resistant. Suitable for all tissue types. Skin friendly yarn with anti-bacterial treatment. High moisture transportation for increased wearing comfort. Perfect fit. Standard and made-to-measure options.



LASTOFA COTTON

Ideal for patients with sensitive skin and for allergy sufferers. Particularly suitable for soft connective tissue. High cotton content of (23%) directly in contact with the skin. Temperature compensating mixture of materials. Perfect fit. Standard and made-to-measure options.



FLAT KNIT

LASTOFA EXTRA

Flat knitted made to measure care for nearly all parts of the body, For the treatment of oedema, burns and scars, High extensive effective pressure. No pressure points due to extra flat seams. High wearing comfort due to 40% cotton content directly on the skin. Skin friendly, thermo and moisture regulating. Optimal fit also in the case of difficult anatomical features. Made-to-measure only.



FOR MORE INFORMATION AND A LOGIN TO THE OFA SHOP CONTACT OAPL

COMPRESSION GARMENTS

SHAPED FULFLEX BELOW KNEE

Shaped tubular bandage that fits below the knee.

Code	Size	Colour	Length	Width
72011B/C	B-C	Beige	1200mm	65-75mm
72012C/D	C-D	Beige	1200mm	70-85mm
72013D/E	D-E	Beige	1200mm	75-85mm
72014E/F	E-F	Beige	120mm	85-100mm



SHAPED FULFLEX ABOVE KNEE

Shaped tubular bandage that fits above the knee.

Code	Size	Colour	Length	Width
72011B/D	B-D	Beige	1200mm	650-750mm
72012C/E	C-E	Beige	1200mm	700-850mm
72013D/F	D-F	Beige	1200mm	750-1000mm
72014E/G	E-G	Beige	120mm	850-1250mm



FULFLEX TUBULAR BANDAGE 1 METRE

A tubular bandage designed to provide compression and support to body segments.

Code	Size	Colour	Width
720101	A	Beige	40mm
720111	B	Beige	65mm
720121	C	Beige	70mm
720131	D	Beige	75mm
720141	E	Beige	85mm
720151	F	Beige	100mm
720161	G	Beige	125mm



FULFLEX TUBULAR BANDAGE 10 METRE

A tubular bandage designed to provide compression and support to body segments.

Code	Size	Colour	Width
72010	A	Beige	40mm
72011	B	Beige	65mm
72012	C	Beige	70mm
72013	D	Beige	75mm
72014	E	Beige	85mm
72015	F	Beige	100mm
72016	G	Beige	125mm
72016H	H	Beige	150mm
72017	J	Beige	175mm
72018	K	Beige	215mm
72019	L	Beige	250mm



BLUE LINE BANDAGE

High quality elastic bandage indicated for the treatment of varicose veins.

Code	Length
72009	3000mm



THE DOFF N' DONNER

The Doff N' Donner takes the work out of applying compression stockings. It conforms to the garment, which increases the life of the garment. Invented with a "push" mentality rather than the "pull" strategy, the Doff N' Donner easily rolls over bandages, wet limbs and even scar management pads. Don & Doff Compression garments of all shapes and compression without frustration.

POLE DOFF N' DONNER

The POLE is the device that "loads" or helps roll the stocking onto the Doff N' Donner. The base of this pole allows for both foot hold as well as suction. We would recommend this item for those that do not have a surface to use a suction cup on.

Please note: The power of the suction cup may leave behind a ring when used on a stained finished surface.



Code

STG-1124
STG-212

Description

Doff n' Donner
Pole Doff n' Donner

VEINOPLUS

To increase blood flow in the legs, reducing venous stasis and minimising the risk of DVT. Conveniently sized stim unit with electrodes included. Suitable for foot and leg oedemas, prevention of DVT, pregnancy, obesity, air and ground travel, prolonged periods of sitting or standing.

Code	Variation	Unit
VEINO	Veinoplus Kit	Each
VEINOE	Replacement Electrodes	P2



JOBST STOCKING DONNER

- Foam grip handles
- Lightweight frame
- Durable baked epoxy finish provides smooth 'no snag' surface

Code

110913



IT STAYS GLUE

- IT STAYS is a roll-on fixative used to hold leg wear in place
- IT STAYS is pliable, moves with your body and can be washed off easily with water
- IT STAYS does not stain

Code	Size
112013	60ml



TERMS & CONDITIONS

ACCOUNT APPLICATION

Trading Name: _____ ACN / ARBN: _____

Legal Name of Entity: _____ ABN: _____

("The Customer")

Applicant's Status: ☐ Company ☐ Sole Trader ☐ Partnership ☐ Trust Government ☐ Other

If individual or partnership, you MUST attach copy of drivers licence(s).

Date of incorporation or date of business started: _____

If change of ownership, date of change: _____

Principal business activity: _____ Expected monthly spend: \$ _____

Registered business address: _____

Delivery address: _____

Contact name: _____ Phone: _____

Mobile: _____ Email: _____

Fax: _____

Account mailing address: _____

Accounts contact person name: _____ Direct Phone: _____

Accounts email : _____ Fax: _____

Preferred method of payment: ☐ Cheque ☐ Credit Card ☐ Bank Transfer (EFT)

Details of Directors/Proprietors:

1) Full Name: _____ DOB: _____

Home Phone: _____ Mobile: _____

Home Address: _____

2) Full Name: _____ DOB: _____

Home Phone: _____ Mobile: _____

Home Address: _____

3) Full Name: _____ DOB: _____
Home Phone: _____ Mobile: _____
Home Address: _____

Trade References

from other commercial suppliers (excluding utilities, telcos etc):

1) Company name: _____ Contact: _____
Phone: _____ Fax: _____
Business Address: _____
2) Company name: _____ Contact: _____
Phone: _____ Fax: _____
Business Address: _____

Banker

1) Bank Name: _____ BSB/Account details _____
12) Bank Name: _____ BSB/Account details _____

I/we acknowledge that: Credit is approved by Orthopaedic Appliances Pty Ltd ("OAPL") in its discretion and that credit may be varied or revoked with written notice; OAPL has informed us and I/we consent to personal and financial information being obtained from and shared with bankers and credit providers in accordance with the Privacy Act 1988 (Cth) and associated laws. I/we consent to the trade and personal references being contacted to gain trading information. I/we the undersigned, accept the Terms & Conditions of Trade which are displayed on OAPL's website (www.oapl.com.au) govern all transactions between me/us and OAPL. I/we agree that any other terms and conditions referred to on any forms that may be used by me/us will apply unless agreed to in writing to the contrary by OAPL. I/we acknowledge that I/we have read the attached Terms and Conditions of Trade and acknowledge that the terms may be varied with written notice. I/we represent to OAPL that this information is true and correct to the best of my/our knowledge.

Acceptance by Authorised person / director / proprietor

Full Name: _____ Position: _____
Signature: _____ Date: _____

GUARANTEE AND INDEMNITY

To: ORTHOPAEDIC APPLIANCES PTY LTD - 29 South Corporate Avenue, Rowville VIC 3178 ("OAPL")

Re: _____ (name of purchasing company [same as on page 1]
ACN _____ [same as on page 1]
Address: _____ (the "Customer").

TERMS & CONDITIONS

1. At the request of the person/s described below (the "Guarantor") OAPL has indicated that it may, at its discretion:
 - a. Sell goods and services under OAPL's standard Terms and Conditions of trade to the Customer.
 - b. Make available certain financial accommodations to the Customer.
2. The Guarantor gives this guarantee and indemnity in consideration of OAPL, and at the request of the Guarantor, providing or continuing to provide financial accommodation to the Customer.
3. The Guarantor unconditionally and irrevocably guarantees to OAPL the due and punctual payment by the Customer to OAPL of all monies (the "Guarantee Monies") at any time actually or contingently owing to OAPL by the Customer jointly and severally on any account, including, without limitation, by way of:
 - a. Monies payable for goods and services purchased by the Customer from OAPL;
 - b. Interest; or
 - c. Costs.
4. If the Customer defaults in payment of the Guarantee Monies, the Guarantor shall pay those monies on demand to, or as directed by, OAPL in order to discharge the Guarantee Monies in full.
5. The demand by OAPL may be sent to the Guarantor at the registered office of the Customer, or on any Guarantor at its address shown below. If posted with postage prepaid, the demand will be taken to be sent three days after posting.
6. The Guarantor unconditionally and irrevocably guarantees to OAPL the due and punctual performance and observance by the Customer of each of OAPL's Terms and Conditions.
7. The Guarantors obligations under this guarantee and indemnity are principal obligations and are not ancillary or collateral to any other right or obligation. The Guarantors obligations may be enforced against the Guarantor without OAPL being required to take action against the Customer or under any security it may hold for the Guarantee Monies
8. The Guarantors obligations under this guarantee and indemnity are a continuing obligation and are not affected by the giving of any time or other indulgence to the Customer by OAPL.
9. This is a continuing guarantee and indemnity. The Guarantors liability is absolute and unconditional and is not affected by anything at all which, but for this provision, might operate to relieve the Guarantor of any obligation in whole or part, including, without limitation:
 - a. Any transaction or arrangement, including but not limited to an arrangement which increases the Guarantee Monies, any extension of time, waiver or release, or any variation in OAPL's Terms and Conditions, which may be made or given between any of OAPL, the Customer, the Guarantor or anyone else;
 - b. The bankruptcy, winding up, liquidation or becoming insolvent under administration of, or appointment of any administrator to the Guarantor, the Customer or anyone else;
 - c. The fact that any other person fails to become bound or ceases to become bound as surety in respect to Guarantee Monies;
 - d. Any failure by OAPL to give notice, or any other omission, delay or mistake on the part of OAPL;
 - e. Any change in the constitution, ownership, directors or management of the Customer, and;
 - f. If the Customer is a trustee of a trust, the Customer acting beyond its powers under the trust.
10. Until the Guarantee Monies have been paid in full, the Guarantor may not exercise any right or subrogation to OAPL and may not exercise any rights of surety in competition with OAPL.
11. The Guarantor unconditionally and irrevocably indemnifies OAPL against any loss or liability which it may suffer because the whole or any part of the Guarantee Monies is not recoverable from the Customer and not recoverable from any Guarantor as a surety, and because of any default by the Customer in the performance and observance by the Customer of each of OAPL's Terms & Conditions.
12. OAPL is at liberty to act as though any one of the Guarantors was the principal debtor and the Guarantors hereby waive all and any of their rights as surety which may at any time be inconsistent with any of these provisions.
13. The Guarantor indemnifies OAPL against, and shall pay OAPL on demand, the amount of all costs, expenses and stamp duty in connections with the enforcement of this guarantee and indemnity.
14. A certificate signed by a director of OAPL or its financial controller setting out the amount owing by the Guarantor pursuant to this guarantee, indemnity and charge at the date specified in the certificate is prima facie evidence in the absence of any manifest error.
15. In the event of any payment of money to OAPL by or on behalf of the Customer or by or on behalf of the Guarantor, being or becoming void, voidable or unenforceable by any statutory provision or otherwise. I/we agree that my/our liability will not have been discharged in respect of such payment and OAPL will have the same rights and remedies against me/us as if that money had not been paid or those transactions had not occurred.

16. The Guarantor is not entitled to set off against or deduct from the Guarantee Monies of any sums owed or claimed to be owed to the Customer by OAPL.
17. If more than one person signs this guarantee and indemnity as Guarantor, they shall be jointly and severally liable.
18. By signing this Guarantee and Indemnity, the Guarantor represents that he/she is not bankrupt and has not been bankrupt and expects that he/she will be able to pay any Guarantee Monies that may be payable under this Guarantee and Indemnity.
19. I give OAPL permission to make inquiries into my financial position and I authorise my bank or other financial institution to provide documents and any information demonstrating my financial position to OAPL upon OAPL's request.
20. I acknowledge that I have read and understood the terms of this Application, Guarantee and Indemnity and Terms and Conditions and that I have received advice from the following Solicitor who explained the contents and meaning of this Application, Guarantee and Indemnity and Terms and Conditions.

Solicitors Firm and Name: _____

on [date]: _____

GUARANTOR (1)

Full Name: _____

Home Address: _____

Signature: _____

Guarantor's Bank name: _____ Guarantor's BSB/ Account: _____

I am the Solicitor for the above mentioned Guarantor. I have been provided with OAPL's Account Application, Guarantee and Indemnity and Terms and Conditions of Trade. In the absence of the Customer and before the Guarantor signed the documents, I explained to the Guarantor the general nature of the documents. In particular I explained that if the Customer defaults in payment to OAPL, the Guarantor will be liable to make good the default which would include all amounts owed, interest and costs. I also explained that giving a guarantee involves substantial risk of losing any security, property or other assets and requires very serious thought. I informed the Guarantor that I was not expressing an opinion or advising on the viability of the transactions with OAPL or the capacity of the Customer or the Guarantor to repay any amount due. At the time of witnessing the Guarantor signing the Guarantee and Indemnity I believe that the Guarantor understood the nature of the Application, Guarantee and Indemnity and the Terms and Conditions and that s/he has signed the agreement of his/her own free will.

SOLICITOR WITNESS

Full Name: _____

Home Address: _____

Signature: _____

GUARANTOR (2)

Full Name: _____

Home Address: _____

Signature: _____

Guarantor's Bank name: _____

Guarantor's BSB/ Account: _____

TERMS & CONDITIONS

I am the Solicitor for the above mentioned Guarantor. I have been provided with OAPL's Account Application, Guarantee and Indemnity and Terms and Conditions of Trade. In the absence of the Customer and before the Guarantor signed the documents, I explained to the Guarantor the general nature of the documents. In particular I explained that if the Customer defaults in payment to OAPL, the Guarantor will be liable to make good the default which would include all amounts owed, interest and costs. I also explained that giving a guarantee involves substantial risk of losing any security, property or other assets and requires very serious thought. I informed the Guarantor that I was not expressing an opinion or advising on the viability of the transactions with OAPL or the capacity of the Customer or the Guarantor to repay any amount due. At the time of witnessing the Guarantor signing the Guarantee and Indemnity I believe that the Guarantor understood the nature of the Application, Guarantee and Indemnity and the Terms and Conditions and that s/he has signed the agreement of his/her own free will.

SOLICITOR WITNESS

Full Name: _____

Home Address: _____

Signature: _____

TERMS & CONDITIONS OF TRADE

1. Definitions

"OAPL" means Orthopaedic Appliances Pty Ltd, ABN 30 006 520 314

"Customer" means the party placing the order with OAPL.

"Order" includes a quotation, order, confirmation of order between OAPL and the Customer which refers to the Goods.

"Goods" means any equipment, services, parts, accessories or materials to be supplied by OAPL.

"Terms" means these Terms and Conditions of Trade, and as amended from time to time.

2. Acceptance of order

2.1 An Order can be made by the Customer in writing or verbally.

2.2 No Order for Goods by the Customer is binding on OAPL until accepted by OAPL. OAPL reserves the right to accept any Order in whole or in part. An individual contract for the supply of Goods is formed by acceptance by OAPL.

2.3 These Terms will override any conditions contained in the Customer's Order. No Order may be cancelled or varied unilaterally by the Customer after acceptance by OAPL.

3. Credit Trading, Financial Information & Price

3.1 The Customer irrevocably authorises OAPL, its servants and agents to make inquiries as they deem necessary to investigate the credit worthiness of the Customer at any time including inquiries upon trade references, Bankers of the Customer, any credit provider to the Customer or credit reporting agency.

3.2 The Customer irrevocably authorises its trade referees, Bankers or other credit provider to disclose to OAPL details of the Customers credit status and financial position as OAPL may require for the purpose of determining the Customers credit worthiness and financial position.

3.3 Credit approval once granted may be withdrawn by OAPL at any time. Where credit approval has been withdrawn, payment for all Goods supplied to the Customer must be paid within 7 days of demand and where no demand is made within 30 days of the withdrawal.

3.4 Where credit has been approved to the Customer, all invoices issued by OAPL are due and payable no later than 30 days after the date of OAPL's invoice or such other date for payment as the Customer and OAPL agree in writing.

3.5 Where credit approval has not been granted, payment of all Goods Ordered by the Customer must be paid before delivery of the Goods to the Customer.

3.6 OAPL reserves the right to change its prices without notice to the Customer.

3.7 Administrative fees may be charged for sales below certain dollar values as determined by OAPL from time to time.

4. Credit card charges and EFT

4.1 If we agree that you can pay by credit card or EFT, then where permitted by law, we reserve the right to levy a fee to cover these facilities.

5. Guarantees under ACL and Exclusions

5.1 The Goods come with guarantees that cannot be excluded under the Australian Consumer Law. The Customer is entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. The Customer is also entitled to have the Goods repaired or replaced if the Goods fail to be of acceptable quality and the failure does not amount to a major failure.

5.2 OAPL's liability in respect of a breach of a guarantee or any warranty made under these Terms for any Goods is limited, to the extent permissible by law and at the option of OAPL to:

- i. The replacement of the Goods or the supply of equivalent goods; or
- ii. The repair of the Goods; or
- iii. Payment of the cost of having the Goods replaced or acquiring equivalent Goods;
- iv. Payment of the cost of having Goods repaired.

5.3 To the extent permitted by law, all other warranties whether implied or otherwise, not set out in these Terms are excluded.

5.4 The Customer indemnifies OAPL and its officers, employees, contractors and agents against any costs, expenses, losses, damages and liability suffered or incurred arising from the Customer's breach of this contract and any negligent or unlawful act or omission of the Customer in connection with the Goods.

6. Warranty

6.1 OAPL provides a warranty against defects within the warranty period after the Goods have been delivered provided that:

- a. The defects have arisen solely from faulty materials or workmanship; and
- b. OAPL is notified of the defect, in writing, within the warranty period; and
- c. The Customer or owner (if the Goods have been transferred) provides proof of purchase or the Goods have been registered with OAPL at the time of purchase; and
- d. The Goods are identified with a label or serial number sufficient to allow OAPL to identify the goods as originating from OAPL, and
- e. The defective Goods are promptly returned to OAPL as required by the Inspection, Claims and Returns clause 7, at no cost to OAPL, and accompanied by a completed Warranty Claim form.

6.2 Excluded from the Warranty are natural wearing parts, wear and tear, damage caused by improper installation, improper loading of weights, misuse, abuse, transportation damage, incorrect or inadequate maintenance, unsuitable location, undue exposure to the elements, unauthorised repairs or alterations carried out by the Customer or third parties, or the use of any spare parts or accessories not manufactured or approved by OAPL.

6.3 Where a replacement is provided, it may not include replacement of the complete product. Items that are beyond economical repair may be given a credit provided that the item has a valid warrantable defect and a replacement is not readily available. OAPL reserves the right to credit the customer in-lieu of replacement; in which case the purchase price of the item will be credited.

7. Inspection, Claims and Returns

7.1 The Customer will inspect and check all Goods received as soon as practicable and must within 48 hours notify OAPL in writing of any shortage in quantity, defect, incorrect specification, damage or Goods not in accordance with the Customer's Order. The Customer is responsible for checking that Goods comply with those Ordered prior to use.

7.2 The Customer agrees that it will not return any of the Goods without first requesting or applying to OAPL so that a Returns Authorisation ("RA") may be issued by OAPL for valid returns.

7.3 OAPL recognises the following categories for valid returns.

- a. Returns accepted within 15 days of invoice date:
 - i. Wrong item delivered;
 - ii. Faulty or damaged product/packaging;
 - iii. Incorrect quantity delivered.
- b. Returns accepted within the warranty period.
- c. The Customer must act promptly to return any Recall items to OAPL.
- d. Returns accepted within 15 days of invoice date:
 - i. Stock item not suitable
 - ii. Order cancelled / no longer required (subject to clause 7.2)
 - iii. Customer ordered incorrectly (subject to clause 7.2)
- e. Returns for Custom Made Goods only if manufactured incorrectly or not in accordance to the Customer's Order. Special orders or non-stock items will not be accepted for return under any circumstances.

7.4 Where an RA has been issued, those Goods to which it related must be returned with ten (10) days of the issue of the RA. The returned Goods must be unmarked, unmodified in a saleable and original condition.

7.5 No credit or replacement will be given unless and until Goods so returned have been verified as being valid for return or are defective and are returned with a valid purchase order/receipt. OAPL reserves the right to recover its costs where returned Goods prove not to be valid for return or are not defective.

7.6 A restocking fee of \$50 per order applies to category Goods referred to in 7.3 d. All freight and handling charges in relation to returning Goods are to be paid by the Customer.

TERMS & CONDITIONS

8. Freight

8.1 Orders of less than \$250 value will incur freight charges. Orders of value greater than \$250 will be shipped freight free.

9. Pricing and Product Descriptions

9.1 Prices are quoted excluding GST unless stated otherwise. Prices are accurate at the time of printing but are subject to change without notice.

9.2 Product descriptions are accurate at the time of printing. Photographs or drawings of products are designed as a reference only. The Customer accepts that the actual product may differ to these.

9.3 OAPL will endeavour to keep Customers informed of the progress of any backorders.

10. Special Orders

10.1 A 40 % deposit is required for all non stock/Special orders prior to the order being placed. This is a non-refundable deposit. All Goods that have been specially ordered for a Customer, may not be returned except where the Customer or the end user is entitled to do so pursuant to the Competition and Consumer Act 2010 ("Cth").

11. Default & Recovery

11.1 In the event of a default by the Customer, OAPL reserves the right to recover all items falling under the terms of the charge in its favour under clause 12.

11.2 The Customer is deemed to be in default immediately upon the happening of any of the following events:

- a. If any payments to OAPL are not made by the due date for payment or payment is not made upon demand by OAPL.
- b. If the Customer ceases to carry on business or stops or suspends payment or states its intention of doing so or is unable to pay its debts as and when they fall due or if any cheque, bill of exchange drawn by the Customer to OAPL is dishonoured.

In the event of a default by the Customer, then without prejudice to any rights which OAPL may have at law or under these Terms:

- a. OAPL or its agents may without notice to the Customer enter the Customer's premises or any premises under the control of the Customer for the purposes of recovering the Goods and the personal property falling under the terms of the charge in its favour under clause 12.
- b. OAPL may recover and resell the Goods and the personal property falling within the terms of the charge in its favour under clause 12.
- c. OAPL will be in no way liable or responsible for any loss or damage to the Goods or for any loss, damage or destruction to the Customer's business howsoever arising from the seizure of the Goods.

12. Security / charge over Goods

12.1 In consideration of the supply of Goods on credit to the Customer from time to time by OAPL and in order to better secure to OAPL all monies payable or to become payable to it pursuant to these Terms, the Customer charges in favour of OAPL all of its right title and interest in any and all Goods supplied by OAPL to it from time to time on credit and which have not been paid for and all proceeds (including but not limited to sale proceeds and an account for such monies) with the amount of the monies owing for such goods.

12.2 The Customer undertakes to keep:

- a. Any Goods supplied by OAPL separate from other goods so that they can be identified as the Goods supplied by OAPL.
- b. Any proceeds in a separate account so that they are identifiable as proceeds arising from a dealing with the charged property referred to herein; and
- c. Records in relation to any receivables account identifying them as proceeds arising under clause 12.1 herein.
- d. Account for the proceeds to OAPL within 7 days of demand.

13. Privacy Policy

13.1 The Customer authorises OAPL to collect the Customers personal information to process the Customers Credit Application. This involves conducting the appropriate checks with your bankers, credit reporting agencies, trade references nominated by you or selected by us, and guarantors of your obligations. For existing trade accounts we may also conduct appropriate checks on an ongoing basis and report bad debtors to the appropriate agency.

13.2 The Customer agrees that OAPL may use the Customers personal information to assist in providing the Goods and Services requested by the Customer, and to improve OAPL's Goods and Services. OAPL or any of its Australian related companies may be in contact with the Customer to inform the Customer about goods, services or promotions which may be of interest to the Customer. We may also share the information OAPL collects with other people or organisations who assist OAPL in providing our goods and services.

13.3 The Customer agrees that OAPL may share the Customers information would like to share your information with these overseas related companies so that they can provide assistance in refining OAPL's Goods and Services, or contact the Customer with special offers. If you would prefer us not to, please write to our Privacy Officer at 93 Webb Street, Fitzroy Vic 3065, facsimile (03 9416 3543, e-mail: privacy@oapl.com.au). You can access your personal information, subject to any exceptions provided by the Privacy Act, by contacting our Privacy Officer.

14. Governing Law

14.1 The Customer agrees that these Terms shall be governed by the laws of Victoria and subject to the non exclusive jurisdiction of the courts of Victoria.

NOTES

Bendigo Clinic

Anne Caudle Centre
Barnard Street
Bendigo VIC 3550
T: 03 5454 8752
F: 03 5454 8756
oapl@bendigohealth.org.au

Box Hill Clinic

Epworth Eastern
Suite 5A, Level 2, Arnold Street
Box Hill VIC 3128
T: 03 9897 1374
F: 03 9783 6944
boxhill@oapl.com.au

Brisbane Clinic

1 / 51 Secam Street
Mansfield QLD 4122
P: 07 3849 8152
F: 07 3849 1354
goodwillortho@oalp.com.au

Bundoora Clinic

Northpark Private Hospital
Greenhills Road
Bundoora VIC 3083
T: 03 9419 2499
F: 03 9416 3543

Brighton Clinic

Suite 4, 214 Bay Street
Brighton VIC 3186
T: 03 9596 6895
F: 03 9824 8205

Canberra Clinic

Mobility Matters
35 Townsville Street
Fyshwick ACT 2609
T: 02 6280 7244
F: 02 6239 1281
oapl@mobilitymatters.com.au

Clayton Clinic

281 Clayton Road
Clayton VIC 3168
T: 03 9562 9422
F: 03 9562 9605
clayton@oapl.com.au

Fitzroy Clinic

93 - 97 Webb Street
Fitzroy VIC 3065
T: 03 9419 2499
F: 03 9416 3543
fitzroy@oapl.com.au

Footscray Clinic

55 Pickett Street
Footscray VIC 3011
T: 03 9687 4589
F: 03 9687 3961
footscray@oapl.com.au

Frankston Clinic

346 Nepean Highway
Frankston VIC 3199
T: 03 9783 3866
F: 03 9783 6944
frankston@oapl.com.au

Perth Clinic

5 Hulme Court
Myaree WA 6154
P: 08 9330 3636
F: 08 9330 3670
info@fos.com.au

Richmond Clinic

Epworth Centre
Suite 5, Level 7, 32 Erin Street
Richmond VIC 3121
T: 03 9421 6226
F: 03 9426 4321
richmond@oapl.com.au

Ringwood Clinic

86 Mt Dandenong Road
Ringwood VIC 3135
T: 03 9879 4299
F: 03 9870 9411
ringwood@oapl.com.au

Sydney Clinic

Macquarie University Clinic 307
Level 3, 2 Technology Place
Macquarie University NSW 2109
P: 02 9812 3854
info@oapl.com.au